



Globe Scientific

95kPA Biohazard Bags

Secure Containment for Safe Air Transport

Globe Scientific 95kPa Biohazard Bags provide secure containment for the safe air transport of liquid specimens. Compliant with 49 CFR 173.27 requirements for transporting dangerous goods by air, they withstand a 95kPa internal pressure differential without rupturing or leaking.

A pressure-sensitive adhesive seal helps contain leaks if the primary specimen container fails, while an integrated outer pouch keeps shipping documents and required paperwork dry and readily accessible.

These bags are ideal for laboratories, healthcare facilities, diagnostic applications, and other settings that require reliable secondary containment for liquid specimens transported by air.



Meets Air Transport Regulations

Compliant with 49 CFR 173.27 for transporting dangerous goods by air. Withstands 95kPa internal pressure differential without rupturing or leaking.



Leak Protection

pressure-sensitive adhesive seal helps contain leaks if the primary specimen container fails.



Keeps Paperwork Dry

Integrated outer pouch keeps shipping documents and required paperwork dry and easily accessible.



Versatile Use

Ideal for laboratories, healthcare facilities, diagnostic applications, and other settings requiring reliable secondary containment for liquid specimens transported by air.

95kPA Biohazard Bags

Secure Containment for Safe Air Transport

Features

- Printed biohazard symbol
- Pressure-sensitive adhesive seal
- Compliant with 49 CFR 173.27 air transport standards
- Manufactured from a proprietary high-barrier multi-layer PE-based composite material
- Integrated external document pouch
- Sealing instructions on front of each bag



Simply fold flap along slit so the star aligns within the box to ensure a proper seal.

Ideal For:



LABORATORIES



HEALTHCARE
FACILITIES



DIAGNOSTIC
APPLICATIONS

Item#	Inner Dimensions	Outer Dimension	Packaging	Unit
4811	6.25" x 8.5"	7" x 10"	50/pack, 14 packs/CS	700
4813	8.25" x 10.5"	9" x 12"	50/pack, 6 packs/CS	300
4814	11" x 13.75"	12" x 15"	50/pack, 5 packs/CS	250
4815	13" x 18"	14" x 19"	50/pack, 2 packs/CS	100

