

Caps for Automated Sample Handling Systems

Seamless cost-effective replacements for OEM caps

Globe Scientific caps for automated systems fit several of the most popular platforms and fit seamlessly into capping processes. These cost-effective alternatives to OEM caps are manufactured to the highest standards to ensure efficient sample processing.

- Seamless cost-effective replacements for OEM caps
- Smooth operation in existing platforms
- Reliable sealing during high-speed transport and agitation
- Made from high-grade Polyethylene (PE) material
- Non-Sterile

13/16mm Cap for **Roche E801**

- Can be used for manual capping of 13/16mm tubes

Item #	Description	Unit
5537	Cap for Roche E80, white, 1000/bag	1000

13mm cap for **IDS / Beckman Track system**

- Can be used for manual capping of 13mm tubes

Item #	Description	Unit
5538	Cap for IDS / Beckman Track system, blue, 500/bag	500

13mm cap for **Beckman DXA5000**

- Can be used for manual capping of 13mm tubes

Item #	Description	Unit
5539	Cap for Beckman DXA5000, clear, 1000/bag	1000



Note: All instruments and names of manufacturers mentioned herein are the trademarks of their respective companies and are used for reference purposes only. Globe Scientific LLC. Expressly disclaims any affiliation to them.