

Calibration Weights

Weigh with confidence when you calibrate your balance with Calibration Weights from Globe Scientific.

Crafted from high-quality, long-lasting stainless steel, our calibration weights come complete with a rugged foam-lined plastic storage case, ensuring accurate calibration for years to come.

With OIML Class E2, F1, and F2 designations, these calibration weights are the perfect companion to any of the Globe Scientific externally calibrating balances.

A Statement of Accuracy comes standard with all Globe Scientific Calibration Weights, or opt for ISO 17025:2017 Calibration performed by an ISO 17025:2017 accredited laboratory.



Calibration Weights

CLASS E2 WEIGHTS – SPECIFICATIONS & ORDERING INFORMATION

Item #	GB-E2-50G	GB-E2-100G	GB-E2-200G
Item # (with Cal. Cert.)	GB-E2-50G-CT	GB-E2-100G-CT	GB-E2-200G-CT
OIML Class	E2	E2	E2
Weight	50 g	100 g	200 g
Tolerance ($\pm$ mg)	± 0.10 mg	± 0.16 mg	± 0.30 mg

CLASS F1 WEIGHTS – SPECIFICATIONS & ORDERING INFORMATION

Item #	GB-F1-100G	GB-F1-200G	GB-F1-500G	GB-F1-1KG	GB-F1-2KG	GB-F1-5KG
Item # (with Cal. Cert.)	GB-F1-100G-CT	GB-F1-200G-CT	GB-F1-500G-CT	GB-F1-1KG-CT	GB-F1-2KG-CT	GB-F1-5KG-CT
OIML Class	F1	F1	F1	F1	F1	F1
Weight	100 g	200 g	500 g	1 kg	2 kg	5 kg
Tolerance ($\pm$ mg)	± 0.5 mg	± 1 mg	± 2.5 mg	± 5 mg	± 10 mg	± 25 mg

CLASS F2 WEIGHTS – SPECIFICATIONS & ORDERING INFORMATION

Item #	GB-F2-500G	GB-F2-1KG	GB-F2-2KG	GB-F2-5KG	GB-F2-10KG	GB-F2-20KG
Item # (with Cal. Cert.)	GB-F2-500G-CT	GB-F2-1KG-CT	GB-F2-2KG-CT	GB-F2-5KG-CT	GB-F2-10KG-CT	GB-F2-20KG-CT
OIML Class	F2	F2	F2	F2	F2	F2
Weight	500 g	1 kg	2 kg	5 kg	10 kg	20 kg
Tolerance ($\pm$ mg)	± 8 mg	± 16 mg	± 30 mg	± 80 mg	± 160 mg	± 300 mg



www.globescientific.com

Customer Service, Sales and Technical Support:

Phone 1 (800) 394-4562 • 1 (201) 599-1400

Fax 1 (201) 599-1406

E-mail mail@globescientific.com

globe
SCIENTIFIC INC.

©2022 Globe Scientific Inc. All rights reserved.

The Globe Scientific logo is a registered trademark of Globe Scientific Inc.