

GENESIS™

BLOOD COLLECTION MIXER

MODEL# CM735A
CAT# 425-CM735A



Our mixer stands alone in the field of blood collection technology – the most accurate, durable and reliable automated mixer.

- Weight accuracy +/- 1 gram
- Monitors blood flow
- Automatic taring
- Reduces QNS rate and over bleeds
- Rechargeable battery
- Perfect for mobile and fixed sites
- Enhanced digital display



GENESIS BPS™

Technology is in our Blood.™

WHEN YOU'RE OUT FOR BLOOD...TAKE US ALONG.

GENESIS™ BLOOD COLLECTION MIXER

MODEL# CM735A

Perfect for fixed and mobile collections
Rechargeable battery or A/C powered



SPECIAL FEATURES

- Rocking action mixes blood with anticoagulant
- Automatically clamps tubing when collected volume reaches pre-set volume
- Enhanced Digital Displays:
 - Pre-set volume by operator
 - Collected volume in mL
 - Total weight in grams
 - Elapsed time of donation
- Audible/Visual Alarms:
 - Slow bleed
 - Low battery
 - End of collection
- Tray:
 - Deeper – Holds all bags & In-line filters
 - Magnetic base ensures stability
 - Leak-proof
 - Tubing guides
- Reduces over and under bleeds (QNS units)
- Automatically tares collection bag & filters
- Rechargeable Ni-Cad battery or 110V/220V A/C power
- Individual recharger included
- Designed for mobile use/transportation case included
- Effective for use in collecting cord blood
- UL/CSA approved
- One year limited warranty



Genesis™ Multi-Charger with Rapid Charge-Discharge Feature

CAT# 425-BC162

- Charge up to 6 batteries simultaneously
- 1 slot for rapid discharge/charge feature
- For maximum battery life, use the rapid discharge feature on each battery once per year.

MODEL# CM735A SPECIFICATIONS

External Dimensions	11.5"L x 8.0"W x 5.5"H
Collection Volume	0–1199 mL
Weight Display	0–1199 grams
Elapsed Time Display	0–19 minutes 59 seconds
Target Volume Range	0–1199 mL Default Value 450 mL or 500 mL with 5 mL increments

Accuracy:

Weight Accuracy	+/- 1 gram
Weight Tolerance	0 grams to + 3 grams

LCD Custom Designed Display

Power Requirements	110V–120V A/C, 220V–240V A/C, 50–60 Hz
--------------------------	---

Battery Pack	Rechargeable Ni-Cad 12V DC
--------------------	----------------------------

Weight:

Mixer with battery	7 lbs
Mixer with case and all accessories	12 lbs
Shipping Weight	17 lbs
Packing	One unit per carton

Accessories Included:

Ni-Cad Battery Pack
Battery Recharger Unit
A/C Power Adapter*
User Manual
Carrying/Transport Case
Mixer Stand

Optional Accessories:

Multi-Battery Charger
Extended Service Agreement

*Suitable for 110V–120V/220V–240V A/C dependent on power supply.

To order, call Genesis BPS at

1-866-712-5663



GENESIS BPS™

Technology is in our Blood.™