

GENESIS™

CENTRIFUGE BALANCE

MODEL# CB220

CAT# 419-CB220



Now it's easy to keep things in balance.

- Electronic differential digital scale
- Displays weight differential of 2 objects
- Balances in seconds
- Simple and easy to operate
- Protects expensive centrifuges from unbalanced loads



GENESIS BPS™

Technology is in our Blood.™

GENESIS™ CENTRIFUGE BALANCE

MODEL# CB220
CAT# 419-CB220



SPECIAL FEATURES

- Electronic differential digital scale
- 2 individual weight sensors and control circuitries
- Independently displays weight differential of 2 objects
- Indicator lamp displays which object weighs less and by how many grams
- Balances in seconds
- Simple and easy to operate
- Protects expensive centrifuges from unbalanced loads

ACCESSORIES INCLUDED WITH PURCHASE OF GENESIS™ CENTRIFUGE BALANCE

- A/C Power adapter
- Spare Fuse 250V/1A
- User Manual

SPECIFICATIONS

Dimensions	13"L x 7"W x 5"H
.....	330L x 180W x 130H mm
Weight measurement	0-2500 grams (range)
Maximum weight difference	2500 grams
Weight sensors	Load cells
Weight display panel	4-digit LED display

INDICATOR LAMPS

Good (if both trays are balanced)
A (left side of balance tray)
B (right side of balance tray)

ACCURACY

Weight calibration tolerance	+/- 1 gram
Decimal truncation error	1 gram
Total accuracy of display	+/- 2 grams
Weight accuracy	4 grams

POWER SOURCE

18V DC (A/C power adapter 100-120/220-240V A/C in)
Power consumption 5 watts

WEIGHT

Centrifuge Balance	7.1 lb
Shipping weight	17.1 lb
Package	one unit per package

To order, call Genesis BPS at

1-866-712-5663



GENESIS BPS™

Technology is in our Blood.™