

# GENESIS™

## CORDLESS AUTOMATIC TUBE STRIPPER

MODEL# BTS-B  
CAT# 49-BTS-B-PRO-MD



### Strip efficiently and extract completely. Anytime. Anywhere.


The GENESIS™ BTS Cordless Automatic Tube Stripper delivers fast and consistent stripping, ensuring complete fluid extraction from tubing. With automatic start-and-stop operation, ergonomic design, and built-in safety protections, the BTS-B streamlines workflow, reduces operator fatigue, and enhances overall processing efficiency.

#### Features:

- Battery operated and rechargeable
- Automatic start/stop function
- Durable, compact benchtop design
- Reduces operator fatigue
- Ensures complete fluid extraction
- Improves workflow efficiency

# GENESIS™ CORDLESS AUTOMATIC TUBE STRIPPER

Strip efficiently. Extract completely. Operate with confidence.



*Tilting rollers help ensure consistent and complete fluid extraction.*

## Also Available:



**Automatic Tube Stripper (BTS)**

**MODEL# BTS**

**CAT# 49-BTS-PRO-MD**

## Ordering Information

|        |                 |
|--------|-----------------|
| Model# | BTS-B           |
| CAT#   | 49-BTS-B-PRO-MD |

## Specifications

|   |   |
|---|---|
| Tube Stripping Speed                    | 27m / min   |
| Dimensions (L x W x H)                  | 10.1" L x 2.1" W x 3.5" H<br>257mm L x 53mm W x 89mm H                      |
| Dimensions (with carrying case)         | 12.6" L x 9.5" W x 2.4" H<br>320mm L x 240mm W x 61mm H                     |
| Weight (without power supply)           | 2.1 lbs / .95 Kg  |
| Weight (with charger and carrying case) | 3.2 lbs / 1.46 Kg   |
| Battery                                 | High Capacity 10.8V 3500mA Battery with Overload and Overcurrent Protection |
| Battery Charger                         | AC/DC IN:100-240V 50-60Hz OUT:14.4V 1A                                      |
| Battery Charging Time                   | Under 6 Hours   |
| Battery Capacity                        | Minimum of 500 Donations (Counted for 1.5m Length Standard Tubing)          |
| Motor                                   | 12V; 25 Watts   |
| Medical Device                          | Class: I  |
| CE Marking                              | Yes   |
| Warranty                                | 1 year  |

**GENESIS** BPS™