

GENESIS™

AUTOMATIC TUBE STRIPPER

MODEL# BTS
CAT# 49-BTS-PRO-MD



Strip efficiently. Extract completely. Operate with confidence.


The GENESIS™ BTS Automatic Tube Stripper delivers fast, hands-free stripping, ensuring consistent and complete fluid extraction from tubing. With automatic start-and-stop operation, ergonomic design, and built-in safety protections, the BTS streamlines workflow, reduces operator fatigue, and enhances overall processing efficiency.

Features:

- Hands-free operation
- Automatic start/stop function
- Durable, compact benchtop design
- Reduces operator fatigue
- Ensures complete fluid extraction
- Improves workflow efficiency

GENESIS™ AUTOMATIC TUBE STRIPPER

Strip efficiently. Extract completely. Operate with confidence.



Tilting rollers help ensure consistent and complete fluid extraction.

Also Available:



Cordless Automatic Tube Stripper (BTS-B)

MODEL# BTS-B

CAT# 49-BTS-B-PRO-MD

Ordering Information

Model#	BTS
CAT#	49-BTS-PRO-MD

Specifications

Tube Stripping Speed	30m / min
Dimensions	7" L x 2.1" W x 3.5" H 179mm L x 53mm W x 89mm H
Dimensions (with carrying case)	12.6" L x 9.5" W x 2.4" H 320mm L x 240mm W x 61mm H
Weight (without power supply)	1.6 lbs / 0.73 Kg
Weight (with power supply & carrying case)	2.8 lbs / 1.27 Kg
Power Supply	AC/DC Adaptor IN:100-240V 50-60Hz OUT:14.4V 1A
Motor	12V; 40 Watts
Medical Device	Class: I
CE Marking	Yes
Warranty	1 year



BTS Table Holder

MODEL# BTS

CAT# 49-BTS-TH-PRO-MD

The BTS Table Holder provides hands-free operation by stabilizing the device for improved control and precision. Its secure design minimizes movement, reduces operator fatigue, and supports a more efficient workflow.

GENESIS BPS™