

GENESIS™

BLOOD COLLECTION MIXER

MODEL# CM735A
CAT# 425-CM735A



Mix with confidence. Collect with precision. The most accurate, durable and reliable automated mixer.

The GENESIS™ Blood Collection Mixer ensures continuous, uniform mixing of blood and anticoagulant for consistent, high-quality collection. Designed for durability and ease of use, it delivers accurate, reliable performance in both fixed and mobile donation settings.

Features:

- Electronically monitors and measures blood flow and volume
- Provides constant mixing of anticoagulant and blood
- Automatically clamps tubing when target volume is achieved
- Automatically tares collection set
- Rechargeable battery and AC power adaptor
- Audible and visual alarms

GENESIS™ BLOOD COLLECTION MIXER

An accurate and reliable automated mixer designed for both fixed site and mobile collections.



LCD Display provides critical collection information.



Automatically clamps tubing at end of collection.



Canvas carrying bag protects all accessories during transport.

Ordering Information Blood Collection Mixer

Model#	CM735A
CAT#	425-CM735A

Specifications

Collection Volume	0-1199 mL
Weight Display	0-1199 Grams
Weight Accuracy	+/- 1 Gram
Weight Tolerance	0 Gram to + 3 Grams
Power Requirements	110V-120V AC 220V-240V AC, 50-60 Hz
Battery Pack	Rechargeable Ni-Cad 12V DC
External Dimension	11.5"L X 8"W X 5.5"H
Weight with Battery	7lbs
Weight with Case & Accessories	12lbs
Shipping Weight	17lbs
Warranty	1 Year (Extended Service Agreements Available)

A-ROM Features

Default Target Volume	500mL
Presets	450mL, 500mL
Alarms	NO 10min and NO 20min Alarms
Clamp	Clamp remains open at 20 minutes

N-ROM Features

Default Target Volume	450mL
Presets	350mL, 400mL, 450mL, 500mL
Alarms	With 10min and 20min Alarms
Clamp	Clamp closes at 20 minutes

Included Components:

- Canvas Carrying Bag
- Magnetic Mounting Tray
- Aluminum Stand
- Rechargeable Ni-Cad Battery
- Individual Battery Charger
- AC Adapter and Power Cord

GENESIS™ Multi-Charger

CAT#: 425-BC162

- Charge up to 6 batteries simultaneously
- 1 Slot for fast charge/discharge feature



GENESIS BPS™