

GENESIS™

CENTRIFUGE BALANCE

MODEL# CB220
CAT# 419-CB220



Dual-scale precision for optimal centrifugation.

The GENESIS™ Centrifuge Balance ensures accurate load matching before centrifugation, reducing vibration and protecting valuable equipment. The intuitive dual digital scale provides instant weight differentials for precise balancing and efficient workflow.

Features:

- Two independent weight sensors with dedicated control circuitry
- Instantly displays weight difference between two containers
- Indicator lamp identifies lighter side and gram variance
- Helps protect centrifuges from unbalanced loads
- Easy to operate, clean and maintain

GENESIS™ CENTRIFUGE BALANCE

Simple and easy to use electronic differential digital scale.



Ordering Information

Model#	CB220
CAT#	419-CB220

Specifications

Dimensions	13" L x 7" W x 5" H 330mm L x 180mm W x 130mm H
Weight Measurement	0-2500g (Range)
Maximum Weight Difference	2500g
Weight Sensors	Load Cells
Weight Display Panel	4-Digit LED Display
Weight Calibration Tolerance	+/- 1g
Decimal Truncation Error	1g
Total Accuracy of Display	+/- 2g
Weight Accuracy	4g
Power Source	18V DC (AC Power Adapter 100-120/220 240V AC in)
Power Consumption	5 Watts
Weight Centrifuge Balance	7.1 lbs / 3.2Kg
Total Shipping Weight	17.1 lbs / 7.8Kg
Warranty	1 year