

# GENESIS™

## RAPID SEAL II™ TUBE SEALER

Bench-top unit with fixed sealing head

MODEL# SE330

CAT# 428-SE330



### High-speed sealing power with proven performance.

The GENESIS™ RapidSeal II™ SE330 is a high-performance RF tube sealer engineered for fast sealing of PVC tubing in blood processing, cell therapy, and biopharmaceutical applications. It produces strong, uniform, easy-to-separate seals in seconds, while its compact, durable benchtop design ensures consistent performance.

#### Features:

- Notched Seal™ for easy separation
- “Intelligent” RF Board – adjusts to all standard PVC blood tubing
- RAPID SEAL II™ seals in 0.5 to 1.2 seconds
- CE and UL compliant

# GENESIS™ RAPID SEAL II™ TUBE SEALER

Bench-top unit with fixed sealing head



## Precision. Reliability. Confidence.

GENESIS™ Tube Sealers provide dependable, high-performance sealing you can trust.

### Ordering Information

Model#	SE330
CAT#	428-SE330

### Specifications

Sealing Head Configuration	Fixed Sealing Head
Dimensions	6.9"W x 6.1"H x 13"D 176mm W x 156mm H x 331mm D
Weight	11.0 lbs / 5.0 Kg
Power Source	AC 100-240V, 50/60 Hz
Max Tubing Diameter	6 mm (1/4")
Warranty	1 Year (Extended Service Agreements Available)

### Included Components:

- Spare Fuses (2)
- Spare Splashguards (5)
- AC Power Cord



*Easy-to-remove head protector simplifies maintenance.*



*Retractable splash guard provides operator safety.*



*Notched Seal™ for easy separation of tubing segments.*



*A complete line of GENESIS™ Tube Sealers is available to meet a wide range of applications and requirements. Please inquire for details.*

**GENESIS** BPS™