

GENESIS™

RAPID SEAL II™ TUBE SEALER

Dual function unit with fixed sealing head
and handheld sealer

MODEL# SE340
CAT# 428-SE340



High-speed sealing power with proven performance.

The GENESIS™ RapidSeal II™ SE340 is a high-performance RF tube sealer engineered for fast sealing of PVC tubing in blood processing, cell therapy, and biopharmaceutical applications. It produces strong, uniform, easy-to-separate seals in seconds using either the fixed head or handheld sealer to ensure consistent performance.

Features:

- Notched Seal™ for easy separation
- "Intelligent" RF Board – adjusts to all standard PVC blood tubing
- RAPID SEAL II™ seals in 0.5 to 1.2 seconds
- 6' standard cable for handheld sealer
- CE and UL compliant

GENESIS™ RAPID SEAL II™ TUBE SEALER

Dual function unit with fixed sealing head and handheld sealer



Precision. Reliability. Confidence.

GENESIS™ Tube Sealers provide dependable, high-performance sealing you can trust.

Ordering Information

Model#	SE340
CAT#	428-SE340

Specifications

Sealing Head Configuration	Fixed Sealing Head and Handheld Sealer
Main Unit Dimensions	6.9"W x 6.1"H x 13"D 176mm W x 156mm H x 331mm D
Handheld Sealing Unit Dimensions	1.1" W x 1.8"H x 8.1"D 28mm W x 44mm H x 204mm D
Main Unit Weight	11.0 lbs / 5.0 Kg
Handheld Weight	0.5 lbs / 0.2 Kg
Power Type	AC
Power Source	AC 100-240V, 50/60 Hz
Max Tubing Diameter	6 mm (1/4")
Handle Cable	6' Standard RF Cable
Warranty	1 Year (Extended Service Agreements Available)

Included Components:

- Handheld Sealer
- 6' Coaxial Cable
- Spare Fuses (2)
- Spare Splashguards (5)
- AC Power Cord



Splash Guard for Hand Sealers

CAT#: 428-SPLGD-HS

- Provides Operator Safety
- Designed for SE340, SE540, and SE640



Easy-to-remove head protector simplifies maintenance.



Retractable splash guard provides operator safety.



Handheld sealer provides operators with control and flexibility.

GENESIS BPS™