

GENESIS™

RAPID SEAL™ SPECIALTY SEALER

Cord Blood and Cellular Therapy
Handheld Sealer

MODEL# SE470
CAT# 428-SE470RS



Advanced sealing for critical cellular applications.

The GENESIS™ RapidSeal™ SE470 is designed for secure, consistent sealing of multi-chamber freezing bags in laboratory and processing environments. It produces strong, uniform, easy-to-separate seals and delivers dependable performance.

Features:

- Notched Seal™ for easy separation
- Specially designed to seal 25 mL Cord Blood Freezing Bag Chambers
- RAPID SEAL™ seals in 1 to 2 seconds
- 6' standard cable for handheld sealer

GENESIS™ RAPID SEAL™ SPECIALTY SEALER

Cord Blood and Cellular Therapy Handheld Sealer



Handheld sealer with elongated jaw, provides operators with control and flexibility.



Sealing Jig

CAT# 428-UCBJIG

- Accurately guides sealing head to seal chambers
- Specially designed to seal 25 mL Cord Blood Freezing Bag Chambers



Precision. Reliability. Confidence.

GENESIS™ Tube Sealers provide dependable, high-performance sealing you can trust.

Ordering Information

Rapid Seal™ Specialty Sealer

Model#	SE470RS
CAT#	428-SE470RS

Ordering Information

Sealing Jig

Model#	UCBJIG
CAT#	428-UCBJIG

Specifications

Sealing Head Configuration	Handheld Sealer
Main Unit Dimensions	4.2"W x 5.9"H x 10.2"D 106mm W x 150mm H x 260mm D
Handheld Sealing Unit Dimensions	1.1" W x 1.8"H x 8.1"D 28mm W x 44mm H x 204mm D
Main Unit Weight	7.7 lbs / 3.5 Kg
Handheld Weight	0.5 lbs / 0.2 Kg
Power Source	AC 100-240V, 50/60 Hz
Coaxial Cable	6' Standard RF Cable
Warranty	1 Year (Extended Service Agreements Available)

Included Components:

- Handheld Sealer
- 6' Coaxial Cable
- Spare Fuses (2)
- Canvas Accessory Bag
- AC Power Cord



A complete line of GENESIS™ Tube Sealers is available to meet a wide range of applications and requirements. Please inquire for details.

GENESIS BPS™