

GENESIS™

RAPID SEAL II™ TUBE SEALER

Battery operated handheld sealer

MODEL# SE640

CAT# 428-SE640



Precision Sealing. Anytime. Anywhere.

The GENESIS™ RapidSeal II™ SE640 is a portable, battery-operated RF tube sealer engineered for fast sealing of PVC tubing in blood processing, cell therapy, and biopharmaceutical applications. It produces strong, uniform, easy-to-separate seals in seconds and supports mobile operations while delivering dependable performance.

Features:

- Battery operated and rechargeable
- 1000+ Seals per full battery charge
- Notched Seal™ for easy separation
- RAPID SEAL II™ seals in 0.5 to 1.2 seconds
- Medical device classification
- 6' standard cable for handheld sealer
- CE and UL compliant

GENESIS™ RAPID SEAL II™ TUBE SEALER

Battery operated handheld sealer.



Precision. Reliability. Confidence.

GENESIS™ Tube Sealers provide dependable, high-performance sealing you can trust.



Handheld sealer provides operators with control and flexibility.



Notched Seal™ for easy separation of tubing segments.



Canvas carrying bag provides easy transport and convenient storage for sealer and all accessories.

Ordering Information

| | |
|--------|-----------|
| Model# | SE640 |
| CAT# | 428-SE640 |

Specifications

| | |
|----------------------------|---|
| Sealing Head Configuration | Handheld Sealer |
| Main Unit Dimensions | 6.9"W x 6.1"H x 13"D 176mm W x 156mm H x 331mm D |
| Handheld Sealer Dimensions | 1.1" W x 1.8"H x 8.1"D 28mm W x 44mm H x 204mm D |
| Main Unit Weight | 11.0 lbs / 5.0 Kg |
| Handheld Weight | 0.5 lbs / 0.2 Kg |
| Power Type | Battery Operated |
| Power Source | 12V / 10AH NiMH Rechargeable Battery |
| Max Tubing Diameter | 6 mm (1/4") |
| Coaxial Cable | 6' Standard RF Cable |
| Warranty | 1 Year (Extended Service Agreements Available) |

Included Components:

- Canvas Carrying Bag
- Handheld Sealer
- 6' Coaxial Cable
- AC Adapter and DC Charger



Splash Guard for Genesis Hand Sealers

CAT:428-SPLGD-HS

- Provides Operator Safety
- Designed for SE340, SE540, and SE640

GENESIS BPS™