

# GENESIS™

## RAPID SEAL II™ TUBE SEALER

Battery operated handheld sealer

MODEL# SE640

CAT# 428-SE640



### Precision Sealing. Anytime. Anywhere.

The GENESIS™ Rapid Seal II™ SE640 is a portable, battery-operated RF tube sealer engineered for fast sealing of PVC tubing in blood processing, cell therapy, and biopharmaceutical applications. It produces strong, uniform, easy-to-separate seals in seconds and supports mobile operations while delivering dependable performance.

#### Features:

- Battery operated and rechargeable
- 1000+ Seals per full battery charge
- Notched Seal™ for easy separation
- RAPID SEAL II™ seals in 0.5 to 1.2 seconds
- Medical device classification
- 6' standard cable for handheld sealer
- CE and UL compliant


# GENESIS™ RAPID SEAL II™ TUBE SEALER

Battery operated handheld sealer.



## Precision. Reliability. Confidence.

GENESIS™ Tube Sealers provide dependable, high-performance sealing you can trust.



*Handheld sealer provides operators with control and flexibility.*



*Notched Seal™ for easy separation of tubing segments.*



*Canvas carrying bag provides easy transport and convenient storage for sealer and all accessories.*

### Ordering Information

Model#	SE640
CAT#	428-SE640

### Specifications

Sealing Head Configuration	Handheld Sealer
Main Unit Dimensions	6.1"W x 1.5"H x 14.4"D 155mm W x 38mm H x 366mm D
Handheld Sealer Dimensions	1.1" W x 1.8"H x 8.1"D 28mm W x 44mm H x 204mm D
Main Unit Weight	11.0 lbs / 5.0 Kg
Handheld Weight	0.5 lbs / 0.2 Kg
Power Type	Battery Operated
Power Source	12V / 10AH NiMH Rechargeable Battery
Max Tubing Diameter	6 mm (1/4")
Coaxial Cable	6' Standard RF Cable
Warranty	1 Year (Extended Service Agreements Available)

### Included Components:

- Canvas Carrying Bag
- Handheld Sealer
- 6' Coaxial Cable
- AC Adapter and DC Charger



### Splash Guard for Genesis Hand Sealers CAT:428-SPLGD-HS

- Provides Operator Safety
- Designed for SE340, SE540, and SE640

**GENESIS** BPS™