

GENESIS™

RAPID SEAL PRO™ TUBE SEALER

Flexible and linkable bench-top sealer with fixed sealing head and optional handheld sealer

MODEL# SE9
CAT# 428-SE9



Durable and reliable high-speed sealing power with proven performance.


The GENESIS™ Rapid Seal Pro™ SE9 is a robust, high-performance tube sealer with a fixed sealing head engineered for fast sealing of PVC tubing in blood processing, cell therapy, biopharmaceutical and commercial applications. It produces strong, uniform, easy-to-separate seals in seconds and provides single-head, linkable multi-head and optional handheld flexibility to maximize efficiency for high volume segmenting.

Features:

- Linkable multi-head configurations
- Stainless steel handheld sealer (optional)
- Notched Seal™ for easy separation
- "Intelligent" RF Board – adjusts to all standard PVC blood tubing
- Rapid Seal Pro™ seals in 0.5 to 1.2 seconds
- 6' standard cable for handheld sealer (optional)
- CE and UL compliant

GENESIS™ RAPID SEAL PRO™ TUBE SEALER

Flexible and linkable bench-top sealer with fixed sealing head and optional handheld sealer



Durable, easy-to-remove metal head protector simplifies maintenance.



Adjustable link bars and thumb screws hold sealing heads firmly in position.



Optional stainless steel handheld sealer provides operators with control and flexibility. (CAT# 428-200-089)

Precision. Reliability. Confidence.

GENESIS™ Tube Sealers provide dependable, high-performance sealing you can trust.



Multi-Head Configurations:

The SE9 Tube Sealer can be used independently or linked as a multi-head sealing unit with flexible head spacing (65/70/75/80/85/90/95/100 mm).



Ordering Information

Model#	SE9
CAT#	428-SE9

Specifications

Sealing Head Configuration	Fixed Sealing Head and Optional Handheld Sealer
Main Unit Dimensions	2.6"W x 6.6"H x 14.3"D 65mm W x 167mm H x 365mm D
Handheld Sealer (Optional) Dimensions	1.0"W x 1.5"H x 7.4"D 25mm W x 39mm H x 188mm D
Main Unit Weight	7.9 lbs / 3.6 Kg
Handheld Sealer Weight (Optional)	0.5 lbs / 0.2 Kg
Power Type	AC
Power Source	AC 100-250V, 50/60 Hz
Max Tubing Diameter	6 mm (1/4")
Coaxial Cable (Optional)	6' Standard RF Cable
Warranty	1 Year (Extended Service Agreements Available)

Included Components:

- Spare Fuses (2)
- Spare Splashguards (5)
- Link Bars (2)
- Thumb Screws (4)
- Signal Cable (1)
- Two-Unit Power Cable
- AC Power Cord



A complete line of GENESIS™ Tube Sealers is available to meet a wide range of applications and requirements. Please inquire for details.

GENESIS BPS™