

Introducing

Diamond™ PEEL & SEAL INOCULATION LOOPS

Precision. Convenience. Control.



A Smarter Way to Access Sterile Loops

Diamond™ Peel & Seal Inoculation Loops simplify microbiology workflows with resealable packaging designed for convenient access and organized storage. When paired with the custom loop stand, the packs remain upright and accessible for efficient, one-handed dispensing under a laminar flow hood.

Ideal Applications

- Microbiology and Clinical laboratories
- Environmental testing
- Food and beverage testing
- Pharmaceutical quality-control laboratories



Diamond™ PEEL & SEAL INOCULATION LOOPS

Designed for Efficient Laboratory Workflows

Resealable Peel & Seal Packaging

The pack can be resealed between uses, helping reduce contamination risk and unnecessary product waste.

Flexible and Rigid Options

Choose flexible loops for greater versatility or rigid loops for enhanced control. Both styles are manufactured from polystyrene (PS).

Available in Two Loop Sizes

Choose from 1 µL and 10 µL loops to meet different sample collection and transfer requirements.

Sterile and Ready to Use


Convenient peel packs contain sterile, ready-to-use loops. The loops are manufactured in an ISO 9001:2015-certified facility and sterilized by electron beam (E-beam) irradiation in accordance with ISO 11137 to a Sterility Assurance Level (SAL) of 10^{-3} .

Available Custom Loop Stand

The proprietary stand securely holds the resealable pack upright, supporting organized storage and convenient one-handed access during repetitive procedures.



Exclusive peel and reseal packs reduce contamination risk and product waste.



Item#	Description	Unit
2881	Loop, Flexible, 1µL, STERILE, Light Green, 40/Peel Pack, 20 Packs/Box	800
2885	Loop, Flexible, 10µL, STERILE, Light Blue, 40/Peel Pack, 20 Packs/Box	800
2891	Loop, Rigid, 1µL, STERILE, Dark Green, 40/Peel Pack, 20 Packs/Box	800
2895	Loop, Rigid, 10µL, STERILE, Dark Blue, 40/Peel Pack, 20 Packs/Box	800
2890	Stand for holding Peel & Seal Inoculation Loop Pack	Each

Packaging

- 40 loops per resealable pack
- 20 packs per box
- 800 loops per box



www.globescientific.com

Phone: 1 (800) 394-4562 • 1 (201) 599-1400 • Fax: 1 (201) 599-1406 • E-mail: mail@globescientific.com

© 2026 Globe Scientific. All rights reserved. The Globe Scientific logo is a registered trademark of Globe Scientific. P&S Loop Flyer 26230