

16 x 101mm Transport Tube with Red Screw Cap, Sterile, Round Bottom, Self-Standing

Item# 6105S-R

Description

Sterile 16 x 101mm polypropylene (PP) transport tube with separate high-density polyethylene (HDPE) red screw cap. These leak-resistant tubes are tested to 95 kPa and E-beam sterilized. The self-standing tubes feature rounded bottoms with skirt.

Features

- Working volume 10mL to 13mL
- Sterilized by Irradiation
- Brim Capacity 15mL
- Self-standing
- Leak resistant
- Non-graduated

Specifications:

- Tubes: USP VI Class approved polypropylene (PP)
- Caps: High-density polyethylene (HDPE)
- Leak resistance tested to 95 kPa
- Cap tightening torque (N·m) 0.50 ±0.05
- Temperature range: -20°C to 120°C
- DNase and RNase free
- Non-Pyrogenic

Quality Control:

- ISO 9001:2015, ISO 13485:2016

Packaging:

Tubes: 200 tubes per individual bag, 5 bags per case

Caps: 1000 caps per bag, 1 bag per case

Products

Item #	Description	Unit
6105S-R	16 x 101mm transport tube, separate red screwcap, STERILE, self-standing, 200 tubes/bag, 1000 caps/case	1000

